

PGE0018 HiBoR Pure

Product Number Packing

PGE0018 : 500G

Product Information

Product Code : PGE0018
Product Name : HiBoR Pure
Shelf Life : 2 years

Technical Specification

Appearance : White crystalline powder
Solubility : Soluble in water
Chloride (Cl) : $\leq 0.001\%$
Calcium (Ca) : $\leq 0.005\%$
Heavy metals (as Pb) : $\leq 0.001\%$
Iron (Fe) : $\leq 0.001\%$
Phosphate (PO₄) : $\leq 0.001\%$
Sulphate (SO₄) : $\leq 0.01\%$
Non volatile with methanol : $\leq 0.05\%$
Insoluble in methanol : $\leq 0.005\%$
Assay : min. 99.50 %

GHS Safety Information

Signal Word : Danger
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 60-61
S-Phrase(s) : 53-45
WGK : 1
RTECS : ED4550000
Storage Temperature(°C) : store below 30°C

Hazard Symbol(s)



Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods