
HilroNTM

PGE0019

Product Identifier

HS Code	:	2833 29 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Bluish-green crystalline powder
Solubility	:	Soluble in water
Chloride (Cl)	:	<= 0.001%
Phosphate (PO ₄)	:	<= 0.001%
Manganese (Mn)	:	<= 0.05%
Copper (Cu)	:	<= 0.005%
Ferric Iron (Fe ⁺³)	:	<= 0.1%
Zinc (Zn)	:	<= 0.005%
Insoluble matter	:	<= 0.01%
Substances not pptd by NH ₄ OH	:	<= 0.05% (as sulphate)
Assay	:	Min. 99.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available