

PGE0020

HiChilCopper

Product Number

Packing

PGE0020 : 100G

Product Information

Product Code : PGE0020
Product Name : HiChilCopper
Shelf Life : 2 years

Technical Specification

Appearance : Blue powder.
Solubility : Soluble in water giving clear blue solution.
Assay : min. 98.00 %

GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 36/37/38
S-Phrase(s) : 26-36
WGK : 3
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods