

PGE0021 HiChilIron

Product Number Packing

PGE0021 : 100G
PGE0021 : 500G

Product Information

Product Code : PGE0021
Product Name : HiChilIron
EC No. : 239-802-2
HS Code : 2922 49 85
Shelf Life : 4 years

Technical Specification

Appearance : Yellowish brown fine crystalline powder
Solubility : Soluble in water
pH (5% in water at 25°C) : 5.5 - 6.5
Lead (Pb) : <= 0.003%
Assay (Fe by Iodometry) : min. 12.00%

GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 36/37/38
S-Phrase(s) : 26-36
WGK : 2
RTECS : AH4900000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

Marine Pollutant : -
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods