

PGE0022**HiMn Pure****Product Number****Packing**

PGE0022	: 500G
PGE0022	: 1KG

Product Information

Product Code	: PGE0022
Product Name	: HiMn Pure
Shelf Life	: 2 years

Technical Specification

Appearance	: Pink crystals or powder
Solubility	: 33.3 mg soluble in 1 mL of water
Chloride (Cl)	: $\leq 0.001\%$
Magnesium (Mg)	: $\leq 0.004\%$
Lead (Pb)	: $\leq 0.001\%$
Iron (Fe)	: $\leq 0.0015\%$
Insoluble matter	: $\leq 0.01\%$
Copper (Cu)	: $\leq 0.001\%$
Sodium (Na)	: $\leq 0.02\%$
Zinc (Zn)	: $\leq 0.01\%$
Potassium (K)	: $\leq 0.01\%$
Nickel (Ni)	: $\leq 0.002\%$
Substances reducing Iodine (I)	: $\leq 0.01\%$
Substances reducing Permanganate	: $\leq 0.001\%$
Oxidising substances	: No reaction
Assay (EDTA Titration)	: min. 98.00 % max. 102.00 %

GHS Safety Information

Hazard Statement(s)	: H373-H411
Precautionary Statement(s)	: P273
Signal Word	: Warning
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s) : 48/20/22-51/53
S-Phrase(s) : 22-61
WGK : 1
RTECS : OP0893500
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

UN No. : 3077
Class : 9
Packaging Group : 3
Marine Pollutant : YES
ADR/RID : 3077 9/PG 3
IMDG : 3077 9/PG 3
IATA : 3077 9/PG 3