



Technical Data

PHM012

Phyto Semiselective Agar Base

Semi-elective medium for the detection of *Xanthomonas campestris* pv. vesicatoria on seeds of pepper and tomato.

Composition **:

Ingredients	Grams/Litre
Soya peptone	2.00
Tryptone	2.00
Glucose anhydrous	1.00
L-glutamine	6.00
L-histidine	1.00
Di-ammonium hydrogen phosphate	0.80
Potassium dihydrogen phosphate	1.00
Magnesium sulphate anhydrous	0.20
Agar	15.00

**Formula adjusted standard to suit the performance parameter

Direction.:

Suspend 29.0 grams in 1000 ml distilled water containing 10 ml Tween 80 . Heat to boiling to dissolve the medium completely. Sterilize the medium by autoclaving at 15 lbs pressure(121°C) for 15 minutes .Cool to 45-50 °C and aseptically add the rehydrated contents of one vial of CCFTNB supplement (PHS011) Mix well and pour into Petri plates .

Principle and Interpretation

Bacterial spot is caused by *Xanthomonas campestris* pv. *vesicatoria* . A semi -selective medium was developed by Sijam (1) wherein *Xanthomonas campestris* pv.vesicatoria was distinguished from its other pathovars by the formation of clear zone around the colony after 2-3 days of incubation.

Soya peptone and tryptone supplies nitrogeous compounds and other nutrients to the organisms. Glucose is the carbohydrate source. L-glutamine and L-histidine supplies amino acid. The phosphates in the medium serves as a buffer. Bacterial and fungal contaminants were inhibited by the antibiotic CCFTNB supplement (PHS011) . 5-fluorouracil in the medium is affective in eliminating fluorescent pseudomonads

PHM012**Phyto Semiselective Agar Base****Quality Control :****Apperance :**

Cream to yellow coloured, homogeneous, free flowing powder.

Gelling :

Firm, comparable with 1.5% Agar gel.

Colour and Clarity of prepared medium :

Yellow coloured, opalescent gel forms in Petri plates

Cultural Response:

Cultural characteristics observed after an incubation at 30-32°C for 5-6 days with added CCFTNB supplement (PHS011) .

Organism (ATCC)	Growth	Colour of the Colony
<i>Xanthomonas campestris</i> pv.vesicatoria	luxuriant	yellow, mucoid, round
<i>Escherichia coli</i> (25922)	inhibited	-
<i>Staphylococcus aureus</i> (25923)	inhibited	-
<i>Pseudomonas aeruginosa</i> (27853)	inhibited	-
<i>Candida albicans</i> (10231)	inhibited	-

References:

1. Sijam, K., Chang, C.J. and Gitaitis, R.D. 1992. A medium for differentiation of tomato and pepper strains of *Xanthomonas campestris* pv. *vesicatoria*. Canad . J. Plant Pathol. 90:208-213

Storage and Shelf-life :

Store below 30°C and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory,diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.