



Technical Data

PHS017

CCFT Supplement

Selective supplement for the detection of *Xanthomonas campestris* pv. *vesicatoria* and *Xanthomonas vesicatoria* on seeds of pepper and tomato.

Formula:

(per vial, sufficient for 1000 ml medium)

Cephalexin monohydrate	65 mg
5-fluorouracil	3 mg
Tobramycin sulphate	0.1 mg
Cycloheximide	100 mg

Directions:

Rehydrate the contents of 1 vial with 2 ml of 0.2N NaOH and 1 ml of methanol. Mix well and aseptically add to sterile molten Phyto Xcv Agar Base (PHM021). Mix well and pour into sterile Petri plates.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.