

PR18 : Optical Sealing Film 96 Well PCR Plate

Product Number : Packing

PR18 : 1X100NO

Product Information

Product Code : PR18

Product Name : Optical Sealing Film 96 Well PCR Plate

Technical Specification

Description : Optical sealing film is designed in such a way to reduce cross contamination between the wells of 96 well plate. The easy-to-use film is applied by peeling away the backing and firmly pressing the sealing film across the top of a microplate. The supreme quality adhesive ensures proper sealing of the plates which helps to reduce sample evaporation during PCR cycles.

Clarity : Clear without any scratches and spots.

Length : 147 mm

Width : 81 mm

Thickness : 0.22 mm

Weight : 2.8 g

DNase/ RNase free : Yes

Applications

- Pressure-sensitive adhesive provides a perfect tight seal over every well of 96 well plate
- Will not interfere with sample reading
- High quality adhesive enables efficient sealing.
- Prevents sample loss and reduces cross-contamination
- Designed to fit 96 and 384 well plates
- These films can be used in various applications including PCR, Real-Time PCR and ELISA.
- Compatible with a wide range of reagents and temperatures

Please refer disclaimer Overleaf.









Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.

Symbols:

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PIPR18

Rev No.: 01

Date of Issue.: 2025-07

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com