

PR27 : Silicone Mat for 96-well plates (PCR)

Product Number : Packing

PR27 : 1X10NO

Product Information

Product Code : PR27

Product Name : Silicone Mat for 96-well plates (PCR)

Technical Specification

Description : Non-sterile, reusable silicone mat with alphanumeric marking

Well Format : 96 well

Capacity : 0.2ml

Dimensions : 115mm x 75mm

Material : Silicone

Autoclavable : Yes

Temperature Range : -80 to 121°C

DNase/ RNase free : Yes

Use : Compatible with 96 Deep well plate and 96 well ELISA plates. These mats reduce evaporation in PCR and storage of samples.

Transport Information

ADR/RID : Not Dangerous Goods


IMDG : Not Dangerous Goods

IATA : Not Dangerous Goods

Please refer disclaimer Overleaf.

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



04/2026

PIPR27_0/0423

PR27-00

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com