

## PCR Blocks, 96-well

**Product code: PR5-1X50NO**

### Description:

- Non-skirted 96-well PCR plates with 0.2 ml volume per well for all standard PCRs and qPCRs. Nuclease-free
- Alpha numeric code for easy sample identification
- This non skirted plate can be sealed using our Sealing Films (PR18/PR20/PR21)

### Quality Control

Each lot of PR5 is assayed for nuclease-free & is functionally tested in DNA/RNA amplification

### Technical Assistance

At HiMedia, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at [mb@himedialabs.com](mailto:mb@himedialabs.com).



Do not use if package is damaged



HiMedia Laboratories Pvt. Limited,  
23 Vadhani Industrial Estate,  
LBS Marg, Mumbai-86, MS, India

PIPR5\_O/0821

PR5-01

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg. office : 23, Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6116 9797 Corporate office : A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)

#### Registered Office :

23, Vadhani Industrial Estate, LBS Marg,  
Mumbai - 400 086, India.  
Tel. : (022) 4017 9797 / 2500 1607  
Fax : (022) 2500 2286

#### Commercial Office

A-516, Swastik Disha Business Park,  
Via Vadhani Indl. Est., LBS Marg,  
Mumbai - 400 086, India  
Tel: 00-91-22-6147 1919  
Fax: 6147 1920, 2500 5764  
Email : [info@himedialabs.com](mailto:info@himedialabs.com)  
Web : [www.himedialabs.com](http://www.himedialabs.com)