

# White Root Culture Modified Medium

With Vitamins and Sucrose  
Without Agar

**Product Code: PT009**

## Product Description :

White Root Culture Modified Medium is a nutrient formulation of macroelements, microelements, vitamins, amino acid and carbohydrate. Potassium nitrate and calcium nitrate serves as a nitrate source and sodium phosphate serves as a phosphate source which helps in the cell division leading to the development of new cells.

Microelements like Boron, Manganese, Molybdenum, Copper, Iron and Zinc enhance metabolism in plants. Boron plays a key role in carbohydrate metabolism. Thiamine, pyridoxine, nicotinic acid act as enzymatic cofactors in universal pathways including glycolysis and TCA cycle along with the primary and secondary metabolism in the plants. Glycine serves as a source of amino acid.

The product is plant tissue culture tested but it is the sole responsibility of the user to ensure the suitability of the medium for individual species.

## Composition :

Ingredients	mg/L
<b>MACROELEMENTS</b>	
Calcium nitrate monohydrate	200.000
Magnesium sulphate	360.000
Potassium nitrate	80.000
Sodium phosphate monobasic	16.500
Sodium sulphate	200.000
<b>MICROELEMENTS</b>	
Boric acid	1.500
Copper sulphate pentahydrate	0.013
Ferrous sulphate heptahydrate	2.500
Manganese sulphate monohydrate	4.500
Molybdic acid (sodium salt)	0.002
Potassium Iodide	0.750
Potassium chloride	65.000
Zinc sulphate heptahydrate	2.800
<b>VITAMINS</b>	
Nicotinic acid (free acid)	0.500
Pyridoxine HCl	0.100
Thiamine hydrochloride	0.100

## AMINO ACID

Glycine 3.000

## CARBOHYDRATE

Sucrose 20000.000

**Total(gms/litre)** 20.9

## Material required but not provided :

- Autoclaved distilled water
- Plant growth regulators
- Gelling agents like Agar (PCT0901) or CleriGel<sup>™</sup>(PCT0903)
- 1N NaOH/HCl

## Precautions :

- Ensure appropriate pH of the medium before addition of gelling agent as acidic pH will lead to decreased gelation resulting in semi solid flowing gel while alkaline pH will lead to formation of hardened gel.
- Use of Distilled water / Tissue culture grade water is recommended for media preparation as tap water or lower grade water may lead to salt precipitation and improper gelation.
- Avoid preparation of concentrated solutions, as it will lead to precipitation of salts.

## Directions :

- Reconstitute medium by adding required quantity of powder in two-third of total volume with constant, gentle stirring till the medium gets completely dissolved.
- Add heat stable supplements prior to autoclaving.
- Make up the final volume with distilled water.
- Adjust the pH of the medium to  $5.75 \pm 0.5$  using 1N NaOH/HCl.
- Add gelling agent and heat the medium to boiling till complete dissolution of gelling agent.
- Sterilize the medium by autoclaving at 15 lbs and 121°C for 15 min.
- Cool the autoclaved medium to about 45°C before adding heat labile supplements.
- Aseptically dispense the desired amount of medium under a laminar airflow unit in sterile culture vessels.

## Quality Control:

### Appearance

White to off-white, homogenous, free flowing powder

### Solubility

20.9 gms/litre soluble in distilled water

### Colour and Clarity

Colourless to light yellow, clear solution

### pH at 25°C

4.40 - 5.40

### Plant Tissue Culture Test

The growth promoting properties of medium is assessed by providing plant cultures with relative humidity of about  $60\% \pm 2\%$ , temperature  $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$  and photoperiod of about 16:8. The plant species showed actively growing callus and shoots with no structural, necrotic and toxic deformity.

## Storage and Shelf Life:

- The plant tissue culture medium powder is extremely hygroscopic and must be stored at  $2-8^{\circ}\text{C}$  in air tight containers.
- Preferably, entire content of each package should be used immediately after opening.
- Use before the expiry date.

Revision : 01 / 2017

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic , research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.