

STERILE DISPOSABLE PETRI PLATES

**PW001/ PW001G/ PW1132/ PW1132G/
PW1132C/ PW1199/ PW1199G/ PW1173**

Transparent, disposable round plates with excellent optical clarity and strength. Petri plate recommended for laboratory use. For invitro use only.

Application : Bioassays, Antimicrobial susceptibility tests , General lab purposes and observation of plant germination.

Product Name	Product Code	Description	Size (mm)	Packing (Nos)
Sterile Disposable Petri Plates	PW001	Polystyrene, optically clear, (ETO Sterilized) transparent disposable Petri Plate	90X15	1X100, 1X600 Individually packed.
Sterile Disposable Petri Plates	PW001G	Polystyrene, optically clear, (Gamma irradiated) transparent disposable Petri Plate	90X15	1X100, 1X600 Individually packed.
Sterile Disposable Petri Plates	PW1132	Polystyrene, optically clear, (ETO Sterilized) transparent disposable Petri Plates	90X15	10X10 60X10 Pack of 10 Plates
Sterile Disposable Petri Plates	PW1132G	Polystyrene, optically clear, transparent disposable Petri Plates (Gamma irradiated)	90X15	10X10 50X10 Pack of 10 Plates
Sterile Disposable Petri Plates	PW1132C	Polystyrene, optically clear, (ETO Sterilized) for easy stacking	90X15	10X10 60X10 Pack of 10 Plates
Sterile Disposable Petri Plates	PW1199	Polystyrene, optically clear, (ETO Sterilized) transparent disposable Petri Plate.	90X15	5X20 30X20 Pack of 20 plates.
Sterile Disposable Petri Plates	PW1199G	Polystyrene, optically clear, (Gamma irradiated) transparent disposable Petri Plates.	90X15	5X20 25X20 Pack of 20 plates.

STERILE DISPOSABLE PETRI PLATES**PW001/ PW001G/ PW1132/ PW1132G/
PW1132C/ PW1199/ PW1199G/ PW1173**

Sterile Disposable Petri Plates Triple Pack.	PW1173	Polystyrene, optically clear, (Gamma irradiated) transparent disposable Petri Plate.	90X15	10X10 50X10 Pack of 10 plates.
--	--------	--	-------	---

Warning :

- Avoid dry heat sterilization.

Disposal : User must ensure safe disposal by autoclaving at a temperature of 121°C (15 lbs) for 15 minutes of used medium (bacteriological materials e.g. culture media with pathogens or spores) in petri plates / unused preparations of this product on completion of work to avoid any contagion.

Product Features :

- Optically Clear for Easy Reading.
- Packaged sterile.
- These Disposable Polystyrene Petri Plates are supplied in a variety of sizes.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane (West) - 400604, Maharashtra, India.
Tel : +91-22-6147 1919 / 6116 9797 / 6903 4800 | Fax : +91-22-6147 1920
Email : info@himedialabs.com | Web : www.himedialabs.com