

## Uricol

PW016 & PW045 .

These containers are available in Sterile and Non-sterile condition They are used for collection of urine sample . Ideal for collecting, storing and transporting urine samples.

**Application :** Medical, Hospital, Infectious Incinerators, Industrial Quality Control, Pathology Laboratories, Academic and Research Institution and Universities.

Product Name	Product Code	Description	Size-Height (mm)	Weight (g)	Capacity (ml)	Materials
Sterile Uricol	PW016	Sterile Uricol with screw capped container for sample collection.	60.10	11.00	50.00	Container- Polystyrene
						Cap - High Density Polyethylene
Non-Sterile Uricol	PW045	Non-Sterile Uricol with screw capped container for sample collection.	60.10	11.00	50.00	Container- Polystyrene
						Cap - High Density Polyethylene

### Product Features :

- Sterile / Non-sterile hard polystyrene container with tight fitting screw cap.
- Quick opening and closing twist cap and leakproof.
- Tight-fitting lids are easy to attach and remove.
- Pack size : 100Nos.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.