

Multipurpose Clinical Sample Collector

PW047 & PW1179.

These are Multipurpose collection container with attached spoon to make collection of samples more convenient . These containers are available in Polystyrene and Polypropylene materials. They are graduated to 100 ml and having spoon attachment. Ideal for collecting, storing and transporting urine samples.

Application : Medical, Hospital, Infectious Incinerators, Industrial Quality Control, Pathology Laboratories, Academic and Research Institution and Universities.

Product Name	Product Code	Description	Size (mm)	Weight (g)	Capacity (ml)	Materials
Multipurpose Clinical Sample Collector	PW047	Autoclavable , graduated and having spoon attachment for sample collection	70 X 54	21.00	100.00	Polypropylene
Sterile Multipurpose Clinical Sample Collector	PW1179	Clear, Sterile , Disposable, graduated and having spoon attachment	71 X 55	24.00	100.00	Polystyrene

Product Features :

- Sample container with spoon
- Material: PP or PS for medical-use
- Quick opening and closing twist cap and leakproof.
- Tight-fitting lids are easy to attach and remove.
- Pack size : 100 Nos.

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.