

## Disposable Petri Plates

**PW001, PW002, PW006, PW007, PW049, PW050, PW054, PW1128, PW1132, PW1132G, PW1134, PW1146, PW1151, PW1155, PW1173, PW1187, PW1190, PW1199, PW1199G, PW1223, PW1226, PW1232 & PW1244**

Transparent, disposable round plates with excellent optical clarity and strength. Petri plate recommended for laboratory use. For invitro use only.

**Application :** Bioassays ,Antimicrobial susceptibility tests , General lab purposes and Observation of plant germination.

Product Name	Product Code	Description	Size (mm)	Packing (Nos)
<b>Sterile Disposable Petri Plate.</b>	<b>PW001</b>	Optically clear, transparent disposable Petri plate	90X15	1X100 , 1X540 & 20X540 Individually packed.
<b>Sterile Disposable Petri Plate.</b>	<b>PW002</b>	Optically clear, transparent disposable Petri plate.	100X15	1X100 & 1X500 Individually packed.
<b>Sterile Disposable Scored Top Petri Plate.</b>	<b>PW006</b>	Optically clear, transparent disposable and useful for culturing work.	90X15	1X100 & 1X500
<b>Sterile Disposable Scored Top Petri Plate.</b>	<b>PW007</b>	Optically clear, transparent disposable and useful for culturing work.	65X12	1X100 & 1X500 Individually packed.
<b>Sterile Disposable Petri Plate.</b>	<b>PW049</b>	Optically clear, transparent disposable and useful for carrying out bioassay.	110X25	1X100 & 1X250
<b>Sterile Disposable Square Petri Plate.</b>	<b>PW050</b>	Optically clear, transparent disposable and useful for carrying out bioassay and antibacterial sensitivity test.	120X120 X15	1X100

<b>Sterile Disposable Petri Plate. (Bi-Petri Plates)</b>	<b>PW054</b>	Divided petri plates. Optically clear, transparent disposable having partition and useful for culturing work.	90X15	1X100 & 1X500
<b>Sterile Disposable Scored Top Petri Plate.</b>	<b>PW1128</b>	Optically clear, transparent disposable petri plates made of polystyrene and useful for culturing work.	75X12	1X100
<b>Sterile Disposable Petri Plate.</b>	<b>PW1132</b>	Optically clear, transparent disposable.	90X15	10 X10, 54 X10 & 1080 X 10 Pack of 10 Plates.
<b>Sterile Disposable Petri Plates.</b>	<b>PW1132G</b>	Optically clear, transparent disposable and gamma irradiated.	90X15	10 X10, 54X10 & 200X10 Pack of 10 plates.
<b>Sterile Disposable Petri Plate.</b>	<b>PW1134</b>	Optically clear, transparent disposable petri-plates.	100X15	10X10 & 50X10 Pack of 10 plates.
<b>Sterile Disposable Petri Plate.</b>	<b>PW1146</b>	Optically clear, transparent disposable and useful for carrying out bioassay and antibacterial sensitivity test.	120X15	1X100 & 1X250
<b>Sterile Disposable Petri Plate.</b>	<b>PW1151</b>	Optically clear, transparent disposable and useful for carrying out bioassay and antibacterial sensitivity test.	110X15	1X100 & 1X250
<b>Sterile Disposable Petri Plate.</b>	<b>PW1155</b>	Optically clear, transparent disposable petri plate.	55X12	5X25 & 10X25. Pack of 25 plates.
<b>Sterile Disposable Petri Plate, Triple Pack.</b>	<b>PW1173</b>	Optically clear, transparent disposable petri plate.	90X15	10X10 & 50X10
<b>Sterile Disposable Scored Top Petri Plate.</b>	<b>PW1187</b>	Optically clear, transparent disposable and useful for culturing work.	55 X 12	1X100 & 1X500
<b>Sterile Disposable Petri Plate.</b>	<b>PW1190</b>	Optically clear, transparent disposable petri plate.	90X25	1X100 & 1X250 Individually

				packed.
<b>Sterile Disposable Petri Plate.</b>	<b>PW1199</b>	Optically clear, transparent disposable petri plate.	90X15	5 X20, 25X20 & 500X20. Pack of 20 plates.
<b>Sterile Disposable Petri Plate.</b>	<b>PW1199G</b>	Optically clear, transparent disposable and gamma irradiated petri plates.	90X15	5 X20 , 25X20 & 500X20 Pack of 20 plates.
<b>Sterile Disposable Petri Plate.</b>	<b>PW1223</b>	Optically clear having 6 vent.	55X14	5X15 & 10X15 Pack of 15 plates.
<b>Sterile Disposable Petri Plate.</b>	<b>PW1226</b>	Optically clear, transparent disposable petri plates. Useful for carrying out bioassay and antibacterial sensitivity test.	150X20	1X100
<b>Sterile Disposable Scored Top Petri Plate.</b>	<b>PW1232</b>	Optically clear, transparent disposable and useful for culturing work.	65X12	10X10 & 50X10 Pack of 10 plates.
<b>Sterile Disposable Petri Plate.(Polystyrene)</b>	<b>PW1244</b>	Optically clear, transparent disposable petri plates made of polystyrene.	110X40	1X100 & 1X250 Individually packed.

**Warning :** Avoid dry heat sterilization.

**Disposal :** User must ensure safe disposal by autoclaving at a temperature of 121°C (15 lbs) for 15 minutes of used medium (bacteriological materials e.g. culture media with pathogens or spores) in petri plates / unused preparations of this product on completion of work to avoid any contagion.

**Product Features :**

- Optically Clear for Easy Reading.
- Packaged sterile.
- These Disposable Polystyrene Petri Dishes are supplied in a variety of sizes.

**Disclaimer :**



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.