

## Disposable Plastic FlexiLoops

**PW012, PW012S, PW013, PW013S, PW014, PW014S, PW051, PW051S, PW052, PW052S, PW053, PW053S, PW1178 & PW1178S.**

Disposable, sterile plastic loops for quantitative procedures such as sampling, serial dilutions and urine counts, as well as for bacterial inoculation. · Molded from inert polymer · Flexible enough to facilitate streaking, yet rigid enough to penetrate the media for subsurface inoculation· Loops available individually wrapped or economical pack ( i.e. 100 loops packed in a steripack container).

**Application :** Microbiology laboratory, Molecular biology laboratory, Plant tissue culture work and various other laboratories.

Product Name	Product Code	Description	Loop Diameter (mm)
Hi-FlexiLoop 2	PW012	Sterilized Flexible loop. Calibrated to 0.005 ml and individually sealed.	2.20
Hi-FlexiLoop 2 Economical Pack	PW012S	Sterilized Flexible loop. Calibrated to 0.005 ml and 100 loops packed in Steripack container.	2.20
Hi-FlexiLoop 4	PW013	Sterilized Flexible loop. Calibrated to 0.01 ml and individually sealed.	4.40
Hi-FlexiLoop 4 Economical Pack	PW013S	Sterilized Flexible loop. Calibrated to 0.01 ml and 100 loops packed in Steripack container.	4.40
Hi-FlexiLoop-S	PW014	Sterilized Flexible loop, Straight for stab inoculation and individually packed.	--
Hi-FlexiLoop-S Economical Pack	PW014S	Sterilized Flexi loop. Straight for stab inoculation and 100 loops packed in Steripack container.	--
StabiFlexiLoop Plus	PW051	Sterilized, Flexible and Calibrated with a streaker / stabber tip . Calibrated to 1µl and individually sealed.	1.25
StabiFlexiLoop Plus Economical Pack	PW051S	Sterilized, Flexible and calibrated with a streaker / stabber tip . Calibrated to 1µl and 100 loops packed in Steripack container.	1.25
StabiFlexiLoop Plus-2	PW052	Sterilized Flexible loop. Calibrated with a streaker / stabber tip. Calibrated to 0.005 ml and individually sealed.	2.20

<b>StabiFlexiLoop-2 Economical Pack</b>	<b>PW052S</b>	Sterilized Flexible loop. Calibrated with a streaker / stabber tip. Calibrated to 0.005 ml and 100 loops packed in Steripack container.	2.20
<b>StabiFlexiLoop Plus-4</b>	<b>PW053</b>	Sterilized Flexible loop. Calibrated with a streaker / stabber tip. Calibrated to 0.01 ml and individually sealed.	4.40
<b>StabiFlexiLoop-3 Economical Pack</b>	<b>PW053S</b>	Sterilized Flexible loop. Calibrated with a streaker / stabber tip. Calibrated to 0.01 ml and 100 loops packed in Steripack container.	4.40
<b>Hi-FlexiLoop 3</b>	<b>PW1178</b>	Sterilized Flexible loop. Calibrated to 0.006 ml and individually sealed.	3.30
<b>Hi-FlexiLoop 3 Economical Pack</b>	<b>PW1178S</b>	Sterilized Flexible loop. Calibrated to 0.006 ml and 100 loops packed in Steripack container.	3.30

#### Product Features:

- They are smooth and flexible to facilitate uniform streaking of culture media.
- Color coded for easy product identification .
- Smooth loop surfaces allows easy streaking.
- Pack Size : 1 X 100 Nos and 5X100 Nos per pack.

#### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.