

HiCaps™

**PW061, PW062, PW063,
PW064 & PW065.**

These caps are designed to fit culture tubes with outer diameters ranging from 12.00mm to 38.00 mm. These Caps are autoclavable and available in various colors.

Application : Microbiology laboratory, Plant tissue culture work and various laboratories.

Product Name	Product Code	Description	Size OD (mm)	Materials
HiCaps-R12	PW061	Assorted colored ,Autoclavable cap	12 .00	Polypropylene
HiCaps-R15	PW062	Assorted colored ,Autoclavable cap	15.00	Polypropylene
HiCaps-R18	PW063	Assorted colored ,Autoclavable cap	18.00	Polypropylene
HiCaps-R25	PW064	Assorted colored ,Autoclavable cap	25.00	Polypropylene
HiCaps-R38	PW065	Assorted colored ,Autoclavable cap	38.00	Polypropylene

Product Features :

- Economical
- Autoclavable .
- Polypropylene are convenient caps for disposable culture tubes.
- Pack size : 100 Nos. and 400Nos.

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.