

Autoclavable Test Tube

PW292 & PW1162.

These tubes are Clear , Transparent, Unbreakable and Rimless. Recommended for the culturing particularly for the media not containing dye.

Application : Used in the laboratories like molecular biology, microbiology, biochemistry ,cell biology and various other laboratories.

Product Name	Product Code	Description	Size (L X D) mm	Capacity (ml)	Material
Autoclavable Test Tube- Rimless	PW292	Clear, Transparent and unbreakable tube for culturing for media not containing dye.	150 X 25	60.00	Polycarbonate
Autoclavable Test Tube - Rimless	PW1162	Clear , Transparent and unbreakable tube for culturing for media not containing dye.	150 X 18.00	30.00	Polycarbonate

Product Features:

- Autoclavable.
- Suitable for media not containing dyes.
- Transparent and unbreakable.
- Made of good quality polycarbonate material.
- Pack size : 100 Nos./Pack

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.