

## Sterile Flocked Nylon Swab With Breakpoint (Scored)

PW1172

Sterile Flocked Nylon Swab with break point (scored) is recommended for nasopharyngeal or urethral specimen collection and viral transport.

### Application:

- Nasopharyngeal or urethral specimen collection
- For Rapid Diagnostics, Medical use
- Bacteriology & Virology

### Size:

Material	Size (Length X Dia)
Shaft	Length : 6.04", Polystyrene
Tip	Diameter :0.135".Length : 0.725".Length of fibre : 0.040" Flocked Nylon fibre

### Qualities :

- Has a breakpoint at which swab can be broken after sampling and put in a tube.
- Nylon fibers attached vertically to shaft, allows for easy collection and release of cells into transport media.
- Small size for pediatric, nasopharyngeal or urethral genital sample collection.

### Product Features:

- Sterile Flocked, Nylon fibre tip
- Flexible Polystyrene shaft
- Pack Size: 1 X 50, 1 X 100 and 1 X 500 Nos.

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.