

Sterile Wooden Stick Without Cotton Bud

PW1176 & PW1177

These multi-purpose sticks are ideal for specimen transfer use, mixing rods and other standard laboratory applications.

Sterile wooden applicator sticks with uniform shape and smooth finish. They can be used for routine spot test sampling and drop transfers, or wherever a handy, disposable mixing aid is required. Also good for irrigating wounds and applying medicines.

Applications: For application specimen transfer, mixing, and more.

Product Name	Product Code	Description	Size (mm)	Packing (Nos).
Sterile Wooden Stick without Cotton bud.	PW1176	Sterile wooden stick without cotton bud. Gamma sterilized and individually packed.	150 X 2.50	1X100 1X500
	PW1177	Sterile wooden stick without cotton bud. Gamma sterilized and pack of 10 stick.	150.0 X 2.50	10 X 10 50 X 10

Product Features

- Multipurpose use.
- Uniform Shape and handy.
- Use for stirring medicines and irrigating wounds
- Also good for applying medicines
- Sterile

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.