

## Sterile Polyester Tipped Swab

PW1180

Sterile Polyester Tipped Swab is recommended for tracheal specimen collection.

### Application:

- Tracheal specimen collection
- Cell Collection for DNA testing
- Gynecological Screening
- Industrial Surface specimen

### Size:

Material	Size ( Length X Dia )
Shaft	6" X 1/10"
Tip	0.687" X 0.203"

### Product Features:

- Sterile regular size tip
- Flexible shaft
- Diagnostic grade polyester
- Pack Size : 1 X 50 and 1X 500 Nos.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.