


PW1211 5 mL Autoclavable Test Tube with Cap Size: 12 X 75mm



Product Number

Packing

PW1211 : 1 X 100 No.

Product Information

Product Code : PW1211
 Product Name : 5 mL Autoclavable Test Tube with Cap Size: 12 X 75 mm
 HS Code : 3926 90 97

Technical Specification

Description : Autoclavable translucent, air tight test tube with cap having conical bottom, self standing and made up of polypropylene material.
 Clarity : Clear without any scratches and fibers.
 Dimension : 12 x 75 mm
 Capacity : 5 mL
 Product alignment : Straight without any bending.
 Temperature tolerance test : Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.

Transport Information

ADR/RID : Not Dangerous Goods
 IMDG : Not Dangerous Goods
 IATA : Not Dangerous Goods