

**PW1247S : Alicote Vial with 'O' Ring (sterile)**

**Product Number : Packing**

PW1247S : 1X1000NO

### Product Information

Product Code : PW1247S

Product Name : Alicote Vial with 'O' Ring (sterile)

### Technical Specification

Description : 1.5ml Screwcap Vial with "O" ring, without hinge, conical bottom, sterile tubes

Clarity : Clear without any particles and fibers.

Color : Clear, natural

Material : Polypropylene

Graduations : No

Nuclease free : Yes

Autoclavable : Yes

Total Height : 49 mm

Body Diameter : 9.8 mm

Inner Diameter (Vial) : 8.2 mm

Over flow capacity (OFC) : 1.9 ml

Weight of the vial : 1.2 g

Temperature Tolerance

Range : -80°C to +121°C

Use : Used to dispense liquid samples

**Please refer disclaimer Overleaf.**

## Transport Information


---

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

## Technical Assistance

---

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at [mb@himedialabs.com](mailto:mb@himedialabs.com).



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,  
Reg. Off: Plot No. C-40, Road No. 21Y,  
MIDC, Wagle Industrial Area, Thane,  
(West) 400604, Maharashtra, INDIA.  
Web: [www.himedialabs.com](http://www.himedialabs.com)



06/2025

PIPW1247S\_1/0324

PW1247S-00

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.  
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)