

PW1292 : Embedding Cassettes, WHITE

Product Number : Packing

PW1292 : 1X1000NO

Product Information

Product Code : PW1292

Product Name : Embedding Cassettes, WHITE

Technical Specification

Description : Tissue Embedding cassettes are made from High density Acetyl Polymer, which is a thermoplastic with high stiffness, low friction and excellent dimensional stability, that keeps specimens safely submerged in solution. These cassettes are totally resistant to most histological solvents.

Color : White

Material : Polyoxymethylene (POM)

Compartment : One

Sterile : No

Autoclavable : Yes

Dimensions : W28mmXH7mmXL42mm

Inner Dimensions (Base) : W26mmXL32mm

Inner Dimensions (Lid) : W24mmXL30mm

Temperature Tolerance

Range : -40°C to +121°C

Use : BIOPSY


Please refer disclaimer Overleaf.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



06/2025

PIPW1292_1/0622

PW1292-01

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg. office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com