

PW1306 : Centrifuge Tube, 15 ml, Orange Cap

Product Number : Packing

PW1306 : 100X1NO
: 4X25NO
: 20X25NO
: 200X25NO

Product Information

Product Code : PW1306
Product Name : Centrifuge Tube, 15 ml, Orange Cap

Technical Specification

Description : Sterile, autoclavable Centrifuge Tube having conical bottom with orange screw cap and printed graduations, large writing area, Nuclease and Pyrogen free

Clarity : Clear without any spots, scratches and fibers.

Color of the Tube : Clear, natural

Color of the Cap : Orange

Material : Polypropylene (Tube)
: Polystyrene (Cap)

Temperature Tolerance

Limit : -20°C to 121°C

Graduations : Yes (0.5 to 15ml)

Parameter : 17000xg max. RCF

Self standing : No

Dimension (ODXL) : 16X120 mm

Cap Diameter : Approx. 30 mm

Please refer disclaimer Overleaf.


Capacity : 15 mL
Product alignment : Straight without any bending

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



09/2025

PIPW1306_2/0922

PW1306-02

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg. office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com