

**PW1309 : Centrifuge Tube, 50 ml, Orange Cap**

**Product Number : Packing**

PW1309 : 100X1NO  
: 4X25NO  
: 20X25NO

**Product Information**

Product Code : PW1309  
Product Name : Centrifuge Tube, 50 ml, Orange Cap

**Technical Specification**

Description : Sterile, Nuclease free, Pyrogen free, Autoclavable Centrifuge Tube, 50ml conical bottom solvent resistant tubes.

Clarity : Clear

Color of the Tube : Natural

Color of the Cap : Orange

Material : Polypropylene (Tube)  
: Polystyrene (Cap)

Temperature Tolerance

Limit : -20°C to 121°C

Graduations : Yes (0.5ml to 3ml, Graduation interval: 0.5ml),  
: Yes (5ml to 50ml, Graduation Interval: 2.5ml)

Autoclavable : Yes

Parameter : 20000xg max. RCF

Self standing : No

Dimension (OD X L) : 30X115 mm

Cap Diameter : Approx. 35mm

**Please refer disclaimer Overleaf.**

Capacity : 50 mL  
 Product alignment : Straight without any bending









**Transport Information**

ADR/RID : Not Dangerous Goods  
 IMDG : Not Dangerous Goods  
 IATA : Not Dangerous Goods

**Technical Assistance**

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at [mb@himedialabs.com](mailto:mb@himedialabs.com).

**Symbols:**

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PIPW1309

Rev No.: 03

Date of Issue.: 2025-08

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg. office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.  
 Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)