

## HiAnticlot

PW131

These vials are for collecting, storing and analyzing blood samples. These vials are sterilized, flat bottom and made of polystyrene material having polypropylene Caps.

Product Name	Product Code	Description	Size (mm)
HiAnticlot Vials Flat Bottom	PW131	Sterile, Polystyrene Vials with Polypropylene Cap.	17 X 54

**Capacity :** Total 6 ml for collection of 3-4 ml blood.

**Application :** For collection of blood sample in Pathology, Hospitals, Medicals and various laboratories & Industries.

### Precautions to be observed:

- Do not autoclave.
- Appropriate for weak acid, weak base and dilute alcohol.

### Product Features:

- Coated with anticoagulant substance.
- 3-4 ml of blood sample can be collected.
- Pack Size : 1 X 500 Nos.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.