

PW1324 : Centrifuge Tube, 50 ml, Orange Cap, self-standing

Product Number : Packing

PW1324 : 100X1NO
: 4X25NO
: 20X25NO

Product Information

Product Code : PW1324
Product Name : Centrifuge Tube, 50 ml, Orange Cap, self-standing

Technical Specification

Description : 50ml Sterile, Nuclease free / Autoclavable, self-standing Centrifuge Tube with conical bottom

Clarity : Clear without any spots, scratches and fibers.

Color of the Tube : Clear, natural

Color of the Cap : Orange

Material : Polypropylene (Tube)
: Polystyrene (Cap)

Graduations : Yes (5 - 50 mL graduation interval: 2.5 mL)

Autoclavable : Yes

Parameter : 5000xg max. RCF

Self- standing : Yes

Dimension (ODXL) : 30X115 mm

Cap Diameter : Approx. 35.2mm

Capacity : 50 mL

Product alignment : Straight without any bending


Please refer disclaimer Overleaf.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



07/2025

PIPW1324_1/0722

PW1324-01

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg. office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com