



Technical Data

PW1327 Screw Tubes, Flat bottom with white PP Cap

Product Number

Packing

PW1327 : 1 X 100 No.

Product Information

Product Code : PW1327
Product Name : Screw Tubes, Flat bottom with white PP Cap

Technical Specification

Description : Clear, tubular, flat bottom glass vial having screw neck made of borosilicate glass, White pp cap with silicone rubber gasket. The entire surface is smooth without presence of any air bubble.

Height of Tube : 99.5 - 100.5 mm

Body Diameter of Tube : 15.70 - 16.70 mm

Inner mouth Diameter of Tube : 9.30 - 9.70 mm

Outer mouth Diameter of Tube : 13.60 - 14.60 mm

Heat resistance test : On exposure to steam sterilization, Tube can withstand heat at 121°C for 30minutes.

Weight with Cap : 15.00 - 16.00 gm

Overflow capacity : 12.0 - 14.0 mL

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Address : C-40, Road No.21/Y, MIDC, Wagle Ind. Area, Thane(W) - 400604, Maharashtra, India. Phone:+91-022-6147 1919, 6903 48 00 Fax:+91-022-6147 1920,6903 48 90 email : info@himedialabs.com. Visit us at our website: www.himedialabs.com