

Tissue Culture Tube With Screw Cap

PW277

Tissue culture tubes with screw cap recommended for tissue culture laboratory use.

Description:

It is a polystyrene, disposable tube. It has conical bottom enabling self standing capability of the tube. The size of the tube is 18mm x 110mm and has a volume capacity of 15.0ml.

Recommended Use:

It can also be used for plant tissue culture work
It can also be used for general purpose laboratory use.

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.