



## HiVeg™ Extract No. 1

RM003V

### Intended use

HiVeg™ Extract No. 1 is prepared under controlled conditions by extracting vegetable proteins. Cultural response is comparable with Meat extract powder (RM003). 0.3 - 0.5% concentration supports the growth of fastidious organisms, when combined with suitable vegetable peptones. It is recommended for general purpose and diagnostic media preparations. It is also suitable for bulk production of antibiotics, enzymes, steroid vaccines and other biological preparations.

### Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Safety guidelines may be referred in individual safety data sheets.

### Limitations

1. It is biological origin product since variation in colour of powder and clarity may observed.
2. Each lot of the product has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's requirement.
3. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium prepared by the product.

### Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature

### Quality Control

- **Appearance** : Light yellow to yellow may have slight green tinge homogeneous free flowing powder characteristic odour but not putrescent
- **Solubility** : Freely soluble in distilled/purified water, insoluble in alcohol.
- **Clarity** : 1% w/v aqueous solution is clear to slight opalescent after autoclaving at 15 lbs pressure (121°C ) for 15 minutes.
- **pH** : pH of 2% w/v aqueous solution at 25°C. 5.50-7.50
- **Microbial Load** :  
Bacterial Count : <= 2000 CFU/gram by plate method, when incubated at 30-35°C for not less than 3 days  
Yeast & mould Count : <= 100 CFU/gram by plate method, when incubated at 20-25°C for not less than 5 days.
- **Test for pathogens** : 1. *Escherichia Coli* - Absent/gram of sample 2. *Salmonella* species - Absent/10 gram of sample 3. *Pseudomonas aeruginosa* - Absent/gram of sample 4. *Staphylococcus aureus* - Absent/gram of sample 5. *Candida albicans* - Absent/gram of sample 6. *Clostridia* - Absent/gram of sample

- **Indole Test** : Tryptophan content: Passes
- **Cultural response** : Cultural response observed after an incubation at 35-37°C for 18-24 hours by preparing Yersinia Isolation HiVeg™ Agar (MV564) using HiVeg™ Extract No.1 as an ingredient.

### Cultural Response

Organism	Growth
<i>Escherichia coli</i> ATCC 25922 (WDCM 00013)	none-poor
<i>Proteus mirabilis</i> ATCC 25933	fair-good
<i>Salmonella enterica</i> subsp. <i>enterica</i> Typhimurium ATCC 14028 (WDCM 00031)	fair-good
<i>Shigella flexneri</i> ATCC 12022 (WDCM 00126)	none-poor
<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i> ATCC 27729	good-luxuriant

### Chemical Analysis :

Total Nitrogen :  $\geq 11.50$  %

Amino Nitrogen :  $\geq 3.50$  %

Sodium chloride :  $\leq 5.00$  %

Loss on drying :  $\leq 7.00$  %

Residue on ignition :  $\leq 15.00$  %

### Storage and Shelf Life

Store between 10 - 30°C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources. Seal the container tightly after use.

### Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Pvt Limited  
C-40,21/Y, MIDC, Wagle Ind Area  
Thane(W)-400604,Maharashtra,India

### Disclaimer :

Revision : 05/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.