

**RM10050 Aniline blue-Orange G Mallory**

**Product Number Packing**

---

RM10050 : 5G

**Product Information**

---

Product Code : RM10050  
Product Name : Aniline blue-Orange G Mallory  
Shelf Life : 4 years

**Technical Specification**

---

Appearance : Orange to dark red powder  
Solubility : 10 ppm in water yields clear solution  
Application : Reagent for a selective collagen staining by aniline blue and for evidencing cytoplasm by Orange G Excellent for modified trichrome staining according to Mallory  
Absorption maxima 1 : 240 - 260 nm (10 ppm in water)  
Absorption maxima 2 : 320 - 340 nm (10 ppm in water)  
Absorption maxima 3 : 480 - 500 nm (10 ppm in water)

**Risk and Safety Information**

---

Storage Temperature(°C) : Store below 30°C

**Transport Information**

---

Marine Pollutant : -  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods