

Chloromethyl isopropyl carbonate

RM10113

Product Identifier

CAS No.	:	35180-01-9
Molecular Formula	:	C ₅ H ₉ ClO ₃
Molecular Weight	:	152.58
HS Code	:	2920 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless liquid.
Solubility	:	33.3 mg soluble in 1 mL of alcohol
Refractive index (n 20/D)	:	1.4110 - 1.4150
Density (at 25°C)	:	1.135 - 1.145 g/mL
Assay (HPLC)	:	min. 99.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Warning
Hazard Statement(s)	:	H315- H317- H411
Precautionary Statement(s)	:	P280
UN No.	:	3082
Class	:	9
Packing Group	:	III
RTECS	:	Not available
WGK	:	3