

**RM10160 N,N-Dimethylacetamide****Product Number Packing**

RM10160 : 500ML

**Product Information**

Product Code : RM10160  
Product Name : N,N-Dimethylacetamide  
Molecular Formula :  $C_4H_9NO$   
Molecular Weight : 87.12  
CAS No. : 127-19-5  
EC No. : 204-826-4  
HS Code : 2924 19 00  
EC Index No. : 616-011-00-4  
Colour Index No. : Laboratory chemicals

**Technical Specification**

Appearance : Colourless liquid.  
Solubility : Soluble in water.  
Density (at 25°C) : ~0.937 g/mL  
Refractive index (n 20/D) : 1.437 - 1.439  
FTIR (Liquid film) : Matches with the standard pattern.  
Assay (GC) : min. 99.00 %

**GHS Safety Information**

Hazard Statement(s) : H312-H332-H360D  
Precautionary Statement(s) : P201-P280-P308 + P313  
Signal Word : Danger  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 61-20/21  
S-Phrase(s) : 53-45  
WGK : 1  
RTECS : AB7700000  
Flash Point(°F) : 158° F  
Flash Point(°C) : 70° C  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**



## Product Information

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.  
email:info@himedialabs.com. Visit us at our website: [www.himedialabs.com](http://www.himedialabs.com)

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein