

RM10179 (-)-Epigallocatechin gallate

Product Number Packing

RM10179 : 10MG

Product Information

Product Code : RM10179
Product Name : (-)-Epigallocatechin gallate
Molecular Formula : $C_{22}H_{18}O_{11}$
Molecular Weight : 458.37
CAS No. : 989-51-5
HS Code : 2932 99 00

Technical Specification

Appearance : White to pink powder.
Solubility : Soluble in water.
FTIR (KBr disc) : Matches with the standard pattern.
Assay (HPLC) : min. 95.00 %

Risk and Safety Information

WGK : 2
RTECS : KB5200000
Storage Temperature(°C) : Store at 2-8°C.

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods