
MUG

RM1019

Product Identifier

CAS No.	:	6160-80-1
EC No.	:	228-186-0
Molecular Formula	:	C ₁₆ H ₁₆ O ₉ .3H ₂ O
Molecular Weight	:	406.33
Synonym	:	4-Methylumbelliferyl-β-D-glucuronide trihydrate
HS Code	:	2932 99 00
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to off-white powder
Solubility	:	33.3 mg soluble in 1 mL of dimethyl sulfoxide
Water (K.F.)	:	12.3 - 14.3%
Assay (HPLC/NaOH T, anhydrous basis)	:	min. 99.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3