

---

## Polyvinylpolypyrrolidone (PVPP)

**RM10324**

### Product Identifier

---

CAS No.	:	9003-39-8
Molecular Formula	:	(C <sub>6</sub> H <sub>9</sub> NO) <sub>n</sub>
Synonym	:	Crospovidone
HS Code	:	3905 99 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White to yellowish-white hygroscopic powder or flakes
Solubility	:	Insoluble in water
Sulphated ash	:	<= 0.10%
Loss on drying (at 100 - 105°C, 2 hr)	:	<= 5.00%
Nitrogen content (dried substance)	:	11.00 - 12.80%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	TR8370000
WGK	:	1