
Sarcosine

RM10341

Product Identifier

CAS No.	:	107-97-1
EC No.	:	203-538-6
Molecular Formula	:	C ₃ H ₇ NO ₂
Molecular Weight	:	89.09
HS Code	:	2922 49 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Melting range	:	205 - 215°C (dec.)
FTIR	:	Matches with the standard pattern.
Assay (NT)	:	97.50 - 102.50%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1