

RM10365

Thioflavine T

Product Number

Packing

RM10365 : 25G

Product Information

Product Code : RM10365
Product Name : Thioflavine T
Synonym : Basic Yellow 1;Thioflavin T
Molecular Formula : C₁₇H₁₉ClN₂S
Molecular Weight : 318.86
CAS No. : 2390-54-7
EC No. : 219-228-9
HS Code : 3212 90 00
Shelf Life : 3 years

Technical Specification

Appearance : Yellow to orange-brown powder
Solubility : 20 ppm in water yields clear solution
Absorption maxima (20 ppm in water): 407 - 417 nm

Risk and Safety Information

WGK : 3
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods