

RM10449

1,2,3,4-Butanetetracarboxylic acid

Product Number

Packing

RM10449 : 100G

Product Information

Product Code : RM10449
Product Name : 1,2,3,4-Butanetetracarboxylic acid
Molecular Formula : C₈H₁₀O₈
Molecular Weight : 234.16
CAS No. : 1703-58-8
EC No. : 216-938-0
HS Code : 2917 19 90
Shelf Life : 4 years

Technical Specification

Appearance : White to off-white crystals or powder
Solubility : 33.3 mg soluble in 1 mL of water
FTIR (Kbr disc) : Matches with the standard pattern
Melting range : 187 - 197°C (dec.)
Assay (NaOH Titration) : 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H302-H319
Precautionary Statement(s) : P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 22-36
S-Phrase(s) : 26
WGK : 3
RTECS : EK6100000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

UN No. : 2503
Class : 8
Packaging Group : 3
Marine Pollutant : No



Product Information

ADR/RID : 2503 8/PG 3
IMDG : 2503 8/PG 3
IATA : 2503 8/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein