

RM10514 Diethylenetriamine**Product Number Packing**

RM10514 : 500ML

Product Information

Product Code : RM10514
Product Name : Diethylenetriamine
Synonym : 2,2'-Diaminodiethylamine
Molecular Formula : C₄H₁₃N₃
Molecular Weight : 103.17
CAS No. : 111-40-0
EC No. : 203-865-4
HS Code : 2921 29 00
EC Index No. : 612-058-00-X
Shelf Life : 4 years

Technical Specification

Appearance : Colourless liquid
Solubility : 1 mL miscible in 1 mL of water
Refractive index (n 20/D) : 1.4826 - 1.4846
FTIR (Liquid film) : Matches with the standard pattern
Density (at 25°C) : ~0.955 g/mL
Assay (GC) : min. 99.00 %

GHS Safety Information

Hazard Statement(s) : H302-H312-H314-H334
Precautionary Statement(s) : P261-P280-P305+P351+P338-P310
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 21/22-34-43
S-Phrase(s) : 26-36/37/39-45
WGK : 2
RTECS : IE1225000
Flash Point(°F) : 215.6 °F
Flash Point(°C) : 102°C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)





Product Information

Transport Information

UN No.	:	2079
Class	:	8
Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	2079 8/PG 2
IMDG	:	2079 8/PG 2
IATA	:	2079 8/PG 2

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein