

RM10543 Finasteride**Product Number Packing**

RM10543	: 100MG
RM10543	: 1G


Product Information

Product Code	: RM10543
Product Name	: Finasteride
Synonym	: N-tert-Butyl-3-oxo-4-aza-5 α -androst-1-en-17 β -carboxamide
Molecular Formula	: C ₂₃ H ₃₆ N ₂ O ₂
Molecular Weight	: 372.54
CAS No.	: 98319-26-7
HS Code	: 2937 29 00
Shelf Life	: 3 years


Technical Specification

Appearance	: White to beige crystals or powder
Solubility	: 33.3 mg soluble in 1 mL of dimethylsulfoxide
Assay (HPLC)	: 98.00 - 102.00%

GHS Safety Information

Hazard Statement(s)	: H302
Signal Word	: Warning
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 22
WGK	: 3
RTECS	: CL5245000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods