
Poloxamer 407

RM10584

Product Identifier

CAS No.	:	9003-11-6
Synonym	:	Pluronic® F-127; Lutrol® F-127
HS Code	:	3402 42 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White wax like granular powder
Solubility	:	10 mg soluble in 1 mL of water
pH (2.5 % in water at 25°C)	:	5.00 - 7.50
Hydroxyl value	:	8 - 11

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	MD0911050
WGK	:	1