

RM10652 Anthraquinone**Product Number Packing**

RM10652 : 500G

Product Information

Product Code : RM10652
Product Name : Anthraquinone
Molecular Formula : $C_{14}H_8O_2$
Molecular Weight : 208.21
CAS No. : 84-65-1
EC No. : 201-549-0
HS Code : 2914 61 00
Shelf Life : 4 years

Technical Specification

Appearance : Yellow to light green to tan crystals or powder
Solubility : 10 mg soluble in 10 mL of chloroform
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 284 - 286°C
Assay (GC/HPLC) : min. 99.00 %

GHS Safety Information

Hazard Statement(s) : H317
Precautionary Statement(s) : P280
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 43
S-Phrase(s) : 36/37
WGK : 3
RTECS : CB4725000
Flash Point(°F) : 365°F
Flash Point(°C) : 185°C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods



Product Information

IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein