

N-(1-Naphthyl)ethylenediamine dihydrochloride, Hi-AR™/ACS

RM1073


Product Identifier

CAS No.	:	1465-25-4
EC No.	:	215-981-2
Molecular Formula	:	C ₁₂ H ₁₄ N ₂ .2HCl
Molecular Weight	:	259.17
HS Code	:	2921 59 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to pink to grey crystals or powder or chunks
Solubility	:	33.3 mg soluble in 1 mL of dilute hydrochloric acid
Sensitivity to diazotised sulphanilamide	:	Passes test
Water (K.F.)	:	<= 5.00%
Assay (AT, on anhydrous basis)	:	98.50 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	KV5330000
WGK	:	3