

## Copper(II) bromide anhydrous

**RM10734**


### Product Identifier

CAS No.	:	7789-45-9
EC No.	:	232-167-2
Molecular Formula	:	CuBr <sub>2</sub>
Molecular Weight	:	223.35
HS Code	:	2827 59 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Black or silver deliquescent crystals or powder
Solubility	:	1000 mg soluble in 1 mL of water
Loss on drying (at 105°C, 2 hr)	:	<= 1.0%
Assay (EDTA Titration, on dry basis)	:	99.00 - 102.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H314- H410
Precautionary Statement(s)	:	P260- P280- P301+P330+P331- P303+P361+P353- P305+P351+P338- P501
UN No.	:	3260
Class	:	8
Packing Group	:	III
RTECS	:	Not available
WGK	:	3