

RM10757 Ethylene glycol diethyl ether**Product Number Packing**

RM10757 : 500ML

Product Information

Product Code : RM10757
Product Name : Ethylene glycol diethyl ether
Molecular Formula : C₆H₁₄O₂
Molecular Weight : 118.17
CAS No. : 629-14-1
EC No. : 211-076-1
HS Code : 2909 19 90
EC Index No. : 603-208-005
Shelf Life : 4 Years

Technical Specification

Appearance : Colourless to pale yellow liquid
Solubility : 33.3 mg soluble in 1 mL of ethanol
FTIR (Liquid film) : Matches with the standard pattern
Refractive index (n 20/D) : 1.391 - 1.395
Density (at 25°C) : 0.837 - 0.847 g/mL
Assay (GC) : min. 98.00%

GHS Safety Information

Hazard Statement(s) : H225-H319-H360Df
Precautionary Statement(s) : P201-P210-P305+P351+P338-P308+P313
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 61-11-19-36-62
S-Phrase(s) : 53-45
WGK : 2
RTECS : KI1225000
Flash Point(°F) : 1.6 °F
Flash Point(°C) : 22 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No



Product Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein