

RM10787 Moxifloxacin hydrochloride**Product Number Packing**

RM10787 : 1G

Product Information

Product Code : RM10787
Product Name : Moxifloxacin hydrochloride
Molecular Formula : $C_{21}H_{24}FN_3O_4.HCl$
Molecular Weight : 437.89
CAS No. : 186826-86-8
HS Code : 2933 49 10
Shelf Life : 3 years

Technical Specification

Appearance : Light yellow to yellow hygroscopic crystals or powder
Solubility : 10 mg soluble in 1 mL of water
FTIR (KBr disc) : Matches with the standard pattern
Water (K.F) : $\leq 4.50\%$
Assay (HPLC on anhydrous basis) : 98.00 - 102.00%

Risk and Safety InformationStorage Temperature($^{\circ}C$) : Store below $30^{\circ}C$ **Transport Information**

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods