

**RM10798 Piperine****Product Number Packing**

RM10798 : 5G

**Product Information**

Product Code : RM10798  
Product Name : Piperine  
Synonym : 1-Piperoylpiperidine  
Molecular Formula : C<sub>17</sub>H<sub>19</sub>NO<sub>3</sub>  
Molecular Weight : 285.34  
CAS No. : 94-62-2  
EC No. : 202-348-0  
HS Code : 2939 99 00

**Technical Specification**

Appearance : White to tan crystals or powder or chunks  
Solubility : 33.3 mg soluble in 1 mL of chloroform  
FTIR (KBr disc) : Matches with the standard pattern  
Melting range : 127 - 135°C  
Assay (HPLC) : min. 98.00%

**GHS Safety Information**

Hazard Statement(s) : H302  
Signal Word : Warning  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22  
S-Phrase(s) : 22-24/25  
WGK : 3  
RTECS : TN2321500  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods