

RM1081 Sodium fluoride, Extra pure**Product Number Packing**

RM2890 : 500G
RM2890 : 5KG

Product Information

Product Code : RM1081
Product Name : Sodium fluoride, Extra pure
Molecular Formula : Na F
Molecular Weight : 41.99
CAS No. : 7681-49-4
EC No. : 231-667-8
HS Code : 2826 19 10
EC Index No. : 009-004-00-7
Shelf Life : 4 years

Technical Specification

Appearance : White crystals or powder
Solubility : 33.3 mg soluble in 1 mL of water
Loss on drying (at 130°C, 3hr.) : <= 0.5%
Chloride (Cl) : <= 0.01%
Iron (Fe) : <= 0.05%
Sulfate (SO4) : <= 0.2%
Assay (NaOH Titration on dry basis) : 98.50 - 102.00%

GHS Safety Information

Hazard Statement(s) : H301-H315-H319
Precautionary Statement(s) : P301+P310-P305+P351+P338
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 25-32-36/38
S-Phrase(s) : 22-36-45
WGK : 1
RTECS : WB0350000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 1690
Class : 6.1



Product Information

Packaging Group	:	3
Marine Pollutant	:	No
ADR/RID	:	1690 6.1/PG 3
IMDG	:	1690 6.1/PG 3
IATA	:	1690 6.1/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein